Sky-Walker

EN	AIR WALKER
PL	BUT MARSZOWY ZE STABILIZACJĄ PNEUMATYCZNĄ
CZ	HLEZENÍ ORTÉZA AIR WALKER
SK	PNEUMATICKÁ ORTÉZA ČLENKA
PT	BOTA IMOBILIZADORA C/AR
EE	JÄIK HÜPPELIIGESE ORTOOS
LV	PĒDAS UN POTĪTES PNEIMATISKAIS FIKSATORS AIRWALKER
LT	ORTOPEDINIS BATAS AIR WALKER
RO	WALKER CU PERNA DE AER
SI	HODILNI ŠKORENJ Z ZRAČNIMI CELICAMI
SE	AIR WALKER BOOT
EL	ΠΟΔΟΚΝΗΜΙΚΌΣ ΝΆΡΘΗΚΑΣ ΜΕ ΑΈΡΑ
HR	ORTOZA AIR WALKER

INSTRUCTIONS FOR USE

Air walker

Product code: **Sky-Walker**

Short air walker

Product code: **Sky-Walker Low**

1. Characteristics

The air walker consists of a splint in the shape of a shoe which is made of PVC, and inner liner of soft material. Owing to the application of modern materials, the orthosis is air permeable and removes excessive heat and moisture perfectly. The air walker is equipped with system of circumferential straps which enable individual adjustment of stabilization force and with profiled anti-sliding sole which facilitates the patient possibly physiological gait. The height of the walker is 17" (Sky-Walker) or 11" (Sky-Walker Low).

2. Application

Air walker is used in post-traumatic and post-operative rehabilitation process in order to limit the range of mobility in ankle joint and to improve functional activities. The orthosis stabilizes the foot and shin in intermediate (neutral) position, in post-traumatic injuries in the area of ankle joint.

Indications for use:

- post-operative treatment
- twist and dislocation of ankle joint
- achilles tendon rupture
- as an alternative for plaster cast immobilization
- fractures in the area of foot, upper and lower ankle joint

3. Contraindications

Do not use in case of:

- skin abrasion
- eczema
- wound and other skin injuries
- after application of warming gel or ointment

Note! In case of varicose veins the usage of the air walker is restricted. Contact with attending physician to determine whether the product may be used together with compression stockings.

4. Notes and recommendations

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

5. Fitting instruction (in sitting position with leg slightly bent in knee joint)

- Unfasten circumferential straps, and take out the inner liner.
- Put the inner liner on the leg. If necessary, use included additional air pads. Fasten the inner liner and put it inside the air walker.
- Put the front shell of the orthosis and fasten circumferential straps by means of Velcro.
- Check the correctness of fitting in standing position. If necessary adjust the force of circumferential straps attachment

Size data (shoe size)

S	36-39 cm
M	40-43 cm
L	44-46,5 cm

6. Maintenance

- Circumferential straps and inner liner can be washed by hand or in washing machine at temperature of 30°C. Do not tumble, chlorinate, iron, bleach, dry cleaning or tumble dry.
- Ensure the velcro is fastened in order not to destroy the material.
- Dry at room temperature far from intensive heat source.
- Clean the PVC splint with damp cloth with soap.
- Dry the splint thoroughly.

7. Storage and transport conditions

Device should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are under a warranty, the terms of which are described in the warranty card received at the time of purchase. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

	Reference number		Manufacturer
	LOT number		Manufacture date
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors.
Learn more at: **www.timago.com**

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